**Legenda dei simboli / Legend of the symbols / Leyenda de los símbolos / Légende des symboles / Legende der Symbole:**

- Leggere le istruzioni per l'uso / Consult instructions for use / Consultar las instrucciones de uso / Consulter le manuel d'utilisation / Die Bedienungsanleitung lesen
- Non riutilizzare / Do not re-use / No reutilizar / Ne pas réutiliser / Nicht wiederverwendbar
- Codice del lotto / Batch code / Código de lote / Code du lot / Chargencode
- Codice del catalogo / Catalog code / Número de catálogo / Référence du catalogue / Produktcode
- Non usare se la confezione è danneggiata / Do not use if package is damaged / No usar si el paquete está dañado / Ne pas utiliser si l'emballage est endommagé / Nicht verwenden, wenn die Verpackung beschädigt ist
- Fabricante / Manufacturer / Fabricante / Fabricant / Hersteller
- Non sterile / Not sterile / No estéril / Nicht steril
- Tenere al riparo dai raggi solari / Keep away from sunlight / Mantener alejado de la luz del sol / Conserver à l'abri du soleil / Vor Sonnenlicht schützen

**IT - CONTENITORI METALLICI, CAPPETTE E CLIP RITENTIVE**

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<b>CLIP GIALLA</b>	±1800gr	<b>CLIP ROSA</b>	±1200gr	<b>CLIP ROSSA</b>	±600gr
--------------------	---------	------------------	---------	-------------------	--------

<b>CLIP BIANCA</b>	<b>STANDARD</b>	<b>CLIP VERDE</b>	<b>SOFT</b>
--------------------	-----------------	-------------------	-------------

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**EN - METAL HOUSINGS, RETENTIVE CAPS AND CLIPS**

**DESCRIPTION OF THE DEVICE AND INDICATIONS OF USE:** the RETENTIVE CAPS, the RETENTIVE CLIPS and the METAL HOUSINGS are spare parts for retentive attachments for the construction of overdentures. **IMPORTANT:** these instructions for use are the most recent available, please read and store them carefully. **IMPORTANT:** These instructions for use are the most recent available; read and store them carefully. **CONTRAINDICATIONS:** at the present state of knowledge there are no contraindications in the use of these products; at the date of emission of these instructions for use no undesirable side effects have been reported. **CAUTION:** the choice of the Medical Device to be used should only be performed by qualified personnel according to the prosthetic design. **MATERIALS:** all the materials used for these products are of quality and meet the international standards of reference. The patient must inform the dentist about potential allergies. If allergic reactions occur search the source and remove it. In case of proven allergy refer to the materials composing the product, do not use them. **WARNINGS AND SINGLE USE:** These products are for **single use only**. Their reuse compromises their mechanical function and can affect the patient's health, as the risk of infection through infectious material transmitted from one patient to another. These products can be used only after washing and sterilization according to the following protocols, which do not compromise their proper functionality. **RHEIN 83 S.R.L.** does not assume any responsibility for any damage resulting from improper use of these products. **TRANSPORT, STORAGE AND EXPIRATION DATE:** do not damage the packaging during transport; store in a dry and clean place in the original boxes away from light and heat sources; the product does not have an expiration date. **WASHING AND STERILIZATION:** the contents of the set are sold as **NOT STERILE**. The metal products must be washed in hot water with a soft brush and non-aggressive and non-corrosive detergent solutions (avoid solutions containing chlorine). The dentist is required to rinse with distilled water, dry, pack and store the prefabricated parts in the original packaging. **STERILIZATION:** the content of the set must be washed in hot water with a soft brush and non-aggressive and non-corrosive detergent solutions (avoid solutions containing chlorine). The dentist is required to rinse with distilled water, dry, pack and store the prefabricated parts in the original packaging. **DISPOSAL:** the disposal of these products must be performed in accordance with the following parameters: steam sterilization with prevacuum cycles at 135° C - 3 minutes, drying time 16 minutes, pressure 2.06 bar. It is recommended to cold disinfect the RETENTIVE CAPS and CLIPS and the PROTECTIVE DISK with benzalkonium chloride solutions for 10 minutes, rinse them with distilled water and dry them before use. The castable parts do not need to be sterilized because of the transformation process (the casting) done by the technician. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided by the manufacturer. The dentist is responsible for the sterilization. **TECHNICAL SUPPORT:** our technical staff is available for technical questions by calling the numbers listed above or by e-mail at: [gianri.storzi@rhein83.it](mailto:gianri.storzi@rhein83.it). For further information on the use of these products, please consult the manual catalogue and the leaflets. **MAINTENANCE AND PERIODIC CARE:** the dentist is responsible for maintaining the correct functionality of the devices with periodic maintenance to ensure the patient's safety. In order to maintain the high-quality standard offered by these products and to avoid the loss of the performances, it is suggested the RETENTIVE CAPS and CLIPS should be replaced periodically according to the instructions provided

<b>COFIAS DORADAS</b> Normo: #500gr. Micro: #450gr.	<b>COFIAS DE PLATA</b> Normo: #350gr. Micro: #300gr.
<b>COFIAS TITAN CAP</b> Normo: #1500gr. Micro: #1300gr.	

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